

Clave 23+ Instruction Manual
Serial numbers starting with 'Y' (4th Gen)
Steam Sterilizer



Distributed by





Instruction Manual
Clave23+

Thank you for choosing our steam sterilizer.

Prior to operating this device, please read the instruction manual carefully and follow all installation instructions.

Regular maintenance:

If this symbol  appears on the screen when the device is powered on or if it shows E88 on the sterilization report, please contact your dealer or local service company for a maintenance service.

This device requires regular maintenance. A  symbol or error code E88 will be displayed after 1 year of use or after a 1000 cycles, whichever comes first. This is programmed into the maintenance screen.

The Clave 23+ is the same as a Clave 23 except it has a built in compressor that speeds up the drying process. It also has a bacteria filter that requires replacement approximately every 6 - 12 months or 1000 cycles.

****USE ONLY DISTILLED WATER IN THIS DEVICE****

Do not re-use the used water from this machine.


Power on the unit with the door open.

This will allow the unit to get an initial zero pressure reading. The machine will not start a cycle and only beep (3 times) when the start button is pressed if it has a non-zero pressure reading on powering up.

General Safety Instructions:

- Make sure that all the installation conditions are fully complied with.
- Ensure that the supply voltage agrees with the supply voltage specified on the rating plate of the sterilizer.
- This appliance must be grounded. Connect only to a properly grounded outlet.
- Do not cover or block any openings on this appliance.
- Use this appliance only for its intended use as described in this manual.
- Do not exceed the maximum weight limit of the loads specified in this manual.
- Do not operate this appliance if it has a damaged power cord or plug or if it is not working properly or if it has been damaged or dropped.
- Never put flammables or explosives products into the sterilizer.
- The sterilizer may not be operated in areas in which gas or any other explosive volatile substance is present.
- Installation and repair work should only be performed by an authorized service technician. Work by unqualified persons could be dangerous and will void any warranty.

Clearing Error Code

Clear Error Code by Pressing and holding the  button. This must be done before running a new cycle.

Document: Version 1.0
Subject to technical changes

Save these instructions

Contents

1. General	4	4 Scope of manual 4 Intended use 4 General safety instructions 4 Standards and directives 5 Symbols
2. Description of the sterilizer	5	5 Sterilizer views 6 Control panel 6 Technical specification 6 Packing content
3. Installation	7	7 General conditions 7 Power supply connection 7 Location requirements
4. Setup	8	8 Fill the distilled water tank 8 Preparation of the sterilization materials 9 Basic set 9 About device 9 Advance set
5. Operation	11	12 Load the sterilization chamber 12 Select the program 12 Start the sterilization program 12 End of cycle 12 Manually interruption of the cycle 13 Test Program 14 Data 14 Save report 14 Printer (optional) 16 Labels (optional)
6. Maintenance	16	16 Clean the distilled water tank 16 Replacement of the bacteriological filter 16 Clean Chamber, trays and tray Rack 17 Door adjustment 17 Replacement of the door seal ring 17 The drain valve
7. Troubleshooting	18	18 Error code description
8. Transportation and storage	19	19 Transportation and storage conditions
9. Safety devices	19	19 Safety device description
11. Appendix	20	20 Water properties/characteristics 21 Diagrams of the sterilization programs 22 Wiring Diagram 23 Hydraulic Diagram

1 General

Scope of Manual

This manual contains information concerning the installation, operation and maintenance of the CLAVE23+ steam sterilizer. To ensure proper performance of the sterilizer, the instructions given in this manual should be thoroughly understood and followed.

Keep the manual near the sterilizer in an accessible location for future reference.

Intended Use

The steam sterilizer described in this manual is intended for the sterilization in all medical, dental, beauty, vet and tattoo fields of the following types of instrument loads: solid, porous, hollow loads type A and hollow loads type B, un-wrapped, single wrapped and double wrapped, and liquid, that are compatible with steam sterilization.

General Safety Instructions

- Read and understand this manual before attempting to install or operate the sterilizer.
- Make sure that all the installation conditions are fully complied with.
- Ensure that the voltage agrees with the supply voltage specified on the rating plate of the sterilizer.
- This appliance must be grounded. Connect only to a properly grounded outlet.
- Do not cover or block any openings on this appliance.
- Use this appliance only for its intended use as described in this manual.
- Do not exceed the maximum weight limit of the loads specified in this manual.
- Do not operate this appliance if it has a damaged cord or plug or if it is not working properly or if it has been damaged or dropped.
- Never put flammables or explosive products into the sterilizer.
- The sterilizer may not be operated in areas in which gas or any other explosive or volatile substance is present.
- Installation and repair work should only be performed by authorized service technicians. Work by unqualified persons could be dangerous and may void the warranty.

Standards and directives

The steam sterilizers were designed and produced in conformity with the following directives and standards:

Directives:

97/23/CE Pressure equipment.

93/42/EEC Medical devices (class II b).

Standards:

EN 13060 Relative to small steam sterilizers.

EN 61010-1 Safety regulations for laboratory devices - Part 1: General regulations.

EN 61010-2-040 Safety regulations specific to sterilizers used in the processing of medical material.

EN 61326-1 Electromagnetic compatibility regulations for laboratory devices.

Symbols

For safe operation, please pay close attention to the alert symbols below which can be found on the sterilizer and throughout this manual.



Important information (Caution)





Hot surface



Ground connection

2 Description of the sterilizer

Sterilizer views











	Date of manufacture
SN	Manufacture number
	Manufacturer Name
	Manufacturer Address
EC REP	EC-Representative Name
	EC-Representative Address

- | | | |
|-------------------------|---|----------------------------|
| 1. Door | 9. Control panel | 16. Circuit breaker |
| 2. Door seal ring | 10. Printer OR Printer Port | 17. Power socket |
| 3. Door lock | 11. USB port | 18. Used water tank vent |
| 4. Door handle | 12. Main switch | 19. Condenser vent |
| 5. Chamber | 13. Distilled water outlet/
Distilled water inlet* | 20. Bacteriological filter |
| 6. Door switch | 14. Used water tank outlet | 21. Distilled water drain* |
| 7. Distilled water tank | 15. Safety valve | 22. Rating plate |
| 8. LCD | | |

*Models equipped with external water filling function.

Control Panel









- 1 Program temperature
- 2 Program
- 3 Holding time
- 4 Pressure
- 5  Distilled water / Fill it flashing.
- 6  Filling the tank
- 7  Bad water quality
- 8  Used water tank full / Drain if flashing
- 9 Printer status
- 10  Door Open
- 11  Door Closed
- 12  Door Locked
- 13 Time
- 14 Current internal temperature
- 15 Up button
- 16 Down button
- 17 Enter button
- 18 Menu button

Technical Specifications

Model	Clave 23+
Chamber (m, diameter/depth)(inches)	φ247 x 450(9.724” x 17.716”)
Overall dimensions (mm,W*H*D)	490 x 455 x 690
Net Weight (kg)	53
Nominal power (VA)	1600
Electrical supply	110V;60 Hz;15A
Sterilization temperatures	121°C/134°C
Capacity of the distilled water tank	2.5 L (Water at level Max.) Approx.
	0.5 L (Water at level Min.) Approx.
Circuit breaker	F15A /400 V
Operation temperature	5°C ~ 40°C
Operation relative humidity	Max. 80%, non-condensing
Max. Noise level	<70 dB
Atmospheric pressure	76 kPa ~ 106 kPa

Packing content

Item	Accessories		Quantity
1	Instrument tray		4
2	Instrument tray rack		1
3	Draining hose		2
4	Door Seal		1
5	Tray handle		1
6	Door adjustment Wrench		1
7	Instruction Manual		1

3 Installation

General conditions

- Position the device on a solid surface with a minimum weight capacity of 60 kgs.
- The sterilizer should be placed on a level worktable.
- Leave at least 10 cm between the device rear part and the wall. The clearance required to open the door is 40 cm.
- Position the sterilizer at such a height as to make it possible for the operator to check the whole sterilization chamber and carry out the normal cleaning operations.
- The room where the device is installed must be sufficiently ventilated.
- Do not install the device near washing basins, taps, etc. where it is likely to be splashed.
- Do not lean on the door when it is opened.
- Do not place trays, papers, fluid containers or other objects on the sterilizer.


Power supply connection

Check the label on the back panel of sterilizer to verify the voltage rating for the unit. Failure to connect the sterilizer to an appropriate power supply could result in damage to the unit and electrical shock to personnel.

Plug the power cord into a properly polarized and grounded receptacle rated for the correct voltage and current. A dedicated circuit that is only used for the sterilizer is recommended. Never connect the device plug to adapters of any type.

Location requirements and dimensions (mm)

To ensure proper air circulation, and to allow access to the reservoir fill port and drain coupling, adhere to the minimum clearance requirements listed below.



4 Setup

Connect the power cord to an outlet of the appropriate voltage. Open the door to remove all of the inner contents for unpacking. Turn on the main power switch on the right side. After switching on, the machine turns on the LCD and shows the door position, water level, working program, date, time, etc.

****THE DOOR MUST BE OPEN WHEN THE POWER IS TURNED ON TO ALLOW CALIBRATION OF THE UNIT.****

Fill the distilled water tank

Manual water filling

When the level of distilled water reaches a minimum level, the distilled water tank icon will flash and beep three times.

Press the button on the tank lid and fully open the door. Fill the tank carefully with distilled water. If the water level exceeds the maximum level, an alarm will sound, and the distilled water tank icon will blink.

Drain the used water tank

Attach the drain hose on the drain port connector located inside the service door. The left port is the drain port for the waste water tank.

Attention: The capacity of the used water tank is approximately 1.5 liters.




Preparation of materials for sterilization


For the most effective sterilization and to preserve instruments, please follow below:

- Clean instruments immediately after use.
- Clean the instruments with an ultrasonic cleaner.
- Residual chemicals left over after the cleaning and disinfecting process may damage and corrode parts of the sterilizer. Always rinse off the instruments using distilled water.
- Follow instrument manufacturer's guidelines and recommendations for handling and cleaning instruments prior to sterilization.
- Check the manufacturer's instructions as to proper procedure for sterilizing each item.
- Arrange the samples of different materials on different trays or with at least 3cm of space between them.
- Clean and dry instruments thoroughly before placing them into a tray.
- Always insert a sterilization paper or cloth between the tray and sample to avoid direct contact.
- Arrange the containers (glasses, cups, test-tubes, etc.) on one side or inverted position, avoiding possible water pooling.
- Don't stack the trays one above the other or put them in direct contact with the walls of the sterilization chamber.
- Always use the instrument tray handle.
- Wrap the samples one by one or, if more tools have to be put in the same bag, verify that these are made of the same material and are not touching one another.
- Don't use metallic clips, pins or similar items, as this jeopardizes the maintenance of the sterilizer.
- Don't overload the trays over the stated limit (see appendix 2).

- Load packaged instruments with the paper side up as the steam and moisture inside the pouches can evaporate easier and thus allow for better drying.


Basic Set


Press the  button to access the main menu. From the main menu, select "Basic Set".






The "Basic Set" menu permits to set the following options:

*Date *Time *Language

Select the "Basic Set" from the main menu by pressing  button.

Select the item by pressing  button. The unit you selected will be highlighted.

Adjust the value by pressing  buttons. Press  button to select the next item.


Press  button to save and exit.

Abbreviation of language options


CHN	Chinese	ENG	English	DEU	German	ESP	Spanish
PL	Polish	FR	French	HUN	Hungary	ROM	Romanian
NL	Dutch	LTU	Lithuanian	LAT	Latvian	CZE	Czech
ITA	Italian	RUS	Russian	PT	Portuguese	HR	Croatian

Note: The Counter (cycle number) cannot be changed by the operator.

About device

Select "About device" from the main menu then press  button.

Press  button to exit.






Advance Set

The "Setup" menu allows you to set the following options:


*Parameter *Unit *Preheat *Expiry date (labels)(Optional) *Water quality (sensor) (optional)

and seeing the information of the: *Last error

Select "Setup" from the main menu by pressing  button.

Input the password one at a time by pressing  and  button to go next.

Password: 1111




Parameter

The “Parameter” menu allows you to set the following options:


*Holding time *Dry time

Select “Parameter” from the menu by pressing button.
Select the program by pressing then press button.
Select the parameter by pressing button. Adjust it by pressing .
Press button to save and exit.



Unit


Select “Unit” from the menu by pressing button.
Select the parameter by pressing button. Adjust it by pressing .
Press button to save and exit.



Preheat


When this mode is activated, the chamber and steam generator start to warm until it reaches the minimum temperature to begin a sterilization cycle. This helps to reduce the total cycle time and improve the drying efficiency. The “Preheat” mode will be deactivated after one hour of inactivity.

Select “Preheat” from the menu by pressing button.
Adjust it by pressing .
Press button to save and exit.




Expiry date (Optional)

To modify the expiration date of the labels, select “Expiry date” from the menu by pressing button. Adjust it by pressing . Press button to save and exit.



Water Quality (optional)


If your sterilizer is equipped with a water quality sensor and you want to deactivate it, select “Water quality” from the menu by pressing button. Adjust it by pressing . Press button to save and exit.



Last Error

In order to help the technical troubleshooting process, the most relevant information corresponding to the last error can be displayed on the screen.

Select “Last error” from the menu by pressing button. Press button to exit.




5 Operation

Prerequisites

Switch On.


Check the status of the icons in the screen, . (Refer to point 2. Description of the sterilizer)

****THE DOOR MUST BE OPEN PRIOR TO TURNING ON THE POWER TO ALLOWING THE UNIT TO CALIBRATE****




Load

Open the door then placed the trays inside the chamber by the tray handle. After the instruments are loaded, you may close the door.




Select the program

Enter to the main menu by pressing button, Select “Program”. Select the program by pressing then press to confirm program, in the screen will appear the information of selected program as the temperature and sterilization time (holding time), also the date, time, current pressure and current temperature.



Start the sterilization program


Press to start the cycle. The stage, conditions and the status of the cycle will appear on the display. The sterilizer will perform the program automatically. (see appendix 2).



Total time or count down until completion ←

End of the cycle

Once the cycle is completed, “End” will appear at the end of the graphic. The printer will print out and the digital report is saved in the USB memory if these are connected.



Caution: Always use the tray handle to load or unload the tray into the sterilizer. Failure to do so can result in burns.

Manual interruption of the cycle

To interrupt a cycle prematurely, hold for 3 seconds.

If the cycle is manually interrupted after it reaches the drying phase, the items inside the sterilizer may be considered sterile and since the cycle has been interrupted during the drying phase the materials and instruments inside the chamber may be wet.

Note: If the cycle is manually interrupted before it reaches the drying phase, the items inside the sterilizer must be considered **not sterile**. N20 will appear on the screen. (see Error code description).



Caution: Caution: Depending which phase of the cycle you interrupt the cycle, steam and water can escape from the sterilization chamber when you open the door.

Data

The internal memory will store the information of the last 9999 cycles.

USB Flash memory (Optional)

A USB drive can be used as a method of storing a report of the cycle. To do so, insert the USB drive into the slot located on the service door of the sterilizer.

The information will automatically output directly to the USB drive after the cycle has completed. The name of the file is determined by the serial number of the machine and the cycle number.

For example:

The serial number is E00001. The cycle number is 0012.

The file name in the USB stick is 01001200.txt.

The first two numbers represent machine number.

The middle four numbers represent cycle number.

The last two numbers represent error code.

E.g. 00:no error; 01: error E01

Printer (Optional)

If installed, you can see the icon in the screen stop flashing.

At the end of each cycle the printer will print out a report of the cycle.


Note: If there is no paper inside the printer, the icon will flash.


Report

Internal Memory

In this menu you can get the information of all the cycles stored in the internal memory of the sterilizer.

Select "Report" from the main menu and press  button, you will see the list of records.

Select the record by pressing  button.

Press  button to print and save the report.


Press  button to exit.

Program Basic Set Report Label

00012 00011 00010 00009

Sample of a printer report

When reading the printed data records, refer to the diagram below:



```
=====
Program: WRAPPED
Temperature: 134C
Pressure: 204.0 kPa
Drying Time: 02Min
Holding Time: 6.0Min
-----
```

```
Time Temp. Pressure
Start 13:21:41 025.0C
T1: 13:27:40 074.8C 051.0kPa
T2: 13:27:57 070.0C 009.7kPa
T3: 13:31:09 088.6C 051.3kPa
T4: 13:31:31 084.0C 009.7kPa
T5: 13:35:42 109.8C 102.7kPa
T6: 13:36:22 098.8C 009.9kPa
TS: 134.8C 222.5kPa
Max. Temperature:135.3C
Min. Temperature:134.2C
Max. Pressure:226.2kPa
Min. Pressure:218.8kPa
T7: 13:44:40 134.3C 218.9kPa
T8: 13:50:39 134.6C 223.1kPa
End 13:53:48 099.0C
-----
```

```
Cycle No.: 00001
Ster. Value: Success
Date: 2018-06-04
SN:E00001
Operator:
v 2N08231C
=====
```

```
=====
Program: WRAPPED
Temperature: 134C
Pressure: 204.0 kPa
Drying Time: 30Min
Holding Time: 6.0Min
-----
```

```
Time Temp. Pressure
Start 15:02:41 023.4C
T1: 15:08:12 072.8C 051.1kPa
T2: 15:08:31 068.6C 009.7kPa
T3: 15:11:30 086.8C 050.5kPa
T4: 15:11:50 083.5C 009.9kPa
T5: 15:15:50 108.2C 101.1kPa
T6: 15:16:30 097.2C 009.7kPa
TS: 000.0C 000.0kPa
Max. Temperature:000.0C
Min. Temperature:000.0C
Max. Pressure:000.0kPa
Min. Pressure:000.0kPa
T7: 15:16:00 000.0C 000.0kPa
T8: 00:00:00 000.0C 000.0kPa
End 15:22:17 116.8C 083.4kPa
-----
```

```
Cycle No.: 00009
Ster. Value: Failure E20
Date: 2018-06-05
SN:E00001
Operator:
v 2N08232C
=====
```

Print labels

Select “Label” from the main menu and press button to enter the menu.
 Select the cycle number by pressing button. Choose the labels quantity by pressing then press button to print.
 Press button to exit.

Program Basic Set Report Label	00012 00011 00010 00009	QT: 01
--	----------------------------------	--------

6 Maintenance

To ensure proper operation and maximum steam sterilizer life, carefully follow all recommendations for periodic maintenance. One of the most important steps you can take to prevent problems with your sterilizer is to use only distilled water.

Frequency	Number of cycles	Maintenance operation
Daily		Clean the external surface
		Drain used water tank
		Clean the door seal
Weekly	25	Clean the chamber, trays and rack
		Clean distilled water tank (drain tank)
Monthly	100	Clean the filter inside the chamber and in the distilled water tank
Every 6 months	500	Replace the bacteria filter
Every year	1000	Replace door seal, tubing, filters, check valves

Clean the distilled water tank

Unplug the power cable, or turn the circuit breaker to the off position at the back of the unit.
 Drain the tank completely using the drain connector at the front of the sterilizer and leave it connected to a tube in an open position.
 Clean the internal surface with a soft sponge and a small soft brush for the areas difficult to reach using distilled water.
 Remove the filter and clean it with a small soft brush and mild soap, rinse it with distilled water, and put it back in to the position.

Replacement of the bacteria filter.

The bacteria filter is at the back of the sterilizer. Unscrew the filter by hand by turning the filter counter-clockwise.
 Replace the bacteria filter with a new one. Screw the new filter by turning the filter clockwise.
 Note: Do not operate sterilizer without filters in place.

Clean chamber, door seal ring, trays, and tray rack.

Remove the trays and tray rack from the chamber. Clean the trays, rack and the inside of the chamber with mild soap.
 Rinse the trays, rack and inside of the chamber with a smooth cloth and distilled water.
 Examine door seal for possible damage. Clean door seal and mating surfaces with a damp cloth.
 Note: Do not use bleaching agents or any abrasive materials or substances in the chamber. Failure to comply may cause damage to the chamber and/or other components.





Caution: To prevent burns, let the unit cool before cleaning gaskets and touching the surface.

Door adjustment

Under normal circumstances, the chamber door does not require adjustments. However, if the seal fails (resulting in steam leaking from the front of the chamber), you may need to adjust the door tension.

Open the door.

Insert the spanner tool in the gap beneath the plastic cover. Use the spanner to grip the adjusting nut. Turn the nut counter clockwise as the figure below. This will tighten the sealing plate. Turn the nut until the sealing plate is tight. If the door knob is too tight, you may also turn the nut clockwise to loosen it.




Caution: Never adjust the chamber door while the door is closed.

Replacement of the door seal ring


Open the chamber door. Remove the door seal ring carefully by hand. Clean the door seal ring carefully with a smooth cloth with distilled water. Moisten the new seal with distilled water. Insert the new seal and tap in sequence as follows:



Press in the top and bottom of the door seal.



Press in the left and right sides of the door seal.




Press the remaining sections of the seal.




Caution: Please ensure the chamber and the door are cold prior to replacing the seal ring.


The drain valves



1. Press the hose onto the drain valve firmly




2. Pull & turn the drain valve clockwise to drain the tank



3. Push & turn the drain valve counter clockwise to close the valve

7 Troubleshooting

Code	Description	Proposed solution
E1	Steam generator temperature sensor error.	Power off & run a new cycle Contact your supplier if error persists.
E2	Inner temperature sensor error.	Power off & run a new cycle Contact your supplier if error persists.
E3	Temperature sensor of the chamber wall error.	Carefully ensure that the chamber wall is heated and contact your supplier.
E5	Fail to release the pressure.	Power off & run a new cycle Contact your supplier if error persists.
E6	Door lock problem during the cycle.	Make sure you had closed the door properly. Check the door switch.
E7	Error between temperature and pressure correlation.	Power off & run a new cycle Contact your supplier if error persists.
E8	Error between temperature and pressure correlation.	Power off & run a new cycle Contact your supplier if error persists.
E9	Failure to hold temperature.	Ensure the distilled tank isn't empty. Check the inner temperature sensor. Check somewhere for leaking.
E10	The door locking system doesn't work.	The electromagnetic locking system doesn't work. The locking system switch doesn't work.
E11	Failure to preheat the steam generator.	Power off & run a new cycle Contact your supplier if error persists.
E12	Failure to preheat the chamber.	Power off & run a new cycle Contact your supplier if error persists.
E16	Pressure error	Replace the air filter Power off & run a new cycle Contact your supplier if error persists.
N20	Program manually interrupted	Reset the error from main screen.
E24	It takes too long time to enter the next status.	Check for leaks. Or contact your supplier if error persists.
E28	The pressure is too high.	Power off and contact your supplier if error persists.
N29	Power Shuts Off During Power Shuts Off During Cycle Cycle	Turn on Power and run new cycle
E34	Pressure is over 50kpa during drying	Power Off and release air from Unit and run another test.

Clear Error Code by Pressing and holding the  button. This must be done before running a new cycle.

8 Transportation and storage

Switch off the sterilizer before transportation or storage.
Pull out the plug. Let the machine cool down.
Drain the distilled water tank and the used water tank.

Conditions for transport and storage
Temperature: -20°C ~ +50°C
Relative humidity: ≤ 85%
Atmospheric pressure: 50kPa~ 106kPa.

9 Safety devices

1. Main breaker: Protection of the instrument against possible failures of the heating elements.
Action: Interruption of the electric power supply.
2. Thermal cutouts on the main transformer winding: protection against possible short circuit and main transformer primary winding overheating
Action: Temporary interruption of winding.
3. Safety valve: Protection against possible sterilization chamber over-pressure.
Action: Release of steam and restoration of the pressure to a safe level.
4. Safety micro-switch for the door status: Comparison for the correct closing position of the door.
Action: Signal of the wrong position of the door.
5. Thermostat on chamber heating elements: Protection for possible overheating of the chamber heating elements.
Action: Interruption of the power supply of the chamber elements.
6. Thermostat on steam generator heating elements: Protection for possible overheating of the steam generator heating elements.
Action: Interruption of the power supply of the steam generator elements.
7. Door safety lock: Protection against accidental opening of the door.
Action: Prevents the accidental opening of the door during the cycle.
8. Self-leveling hydraulic system: Hydraulic system for the natural pressure leveling in case of manual cycle interruption, alarm or blackout.
Action: Automatic restoration of the atmospheric pressure inside the chamber. (valve opens to release internal pressure)

Appendix 1

Water properties / Characteristics

Description	Feed water	Condensate
Evaporate residue	$\leq 10\text{mg/ I}$	$\leq 1.0\text{mg/kg}$
Silicon oxide sio2	$\leq 1\text{mg/ I}$	$\leq 1.0\text{mg/kg}$
Iron	$\leq 0.2\text{mg/ I}$	$\leq 0.1\text{mg/kg}$
Cadmium	$\leq 0.005\text{mg/ I}$	$\leq 0.05\text{mg/kg}$
Lead	$\leq 0.05\text{mg/ I}$	$\leq 0.1\text{mg/kg}$
Rest of heavy metals	$\leq 0.1\text{mg/ I}$	$\leq 0.1\text{mg/kg}$
Chloride	$\leq 2\text{mg/ I}$	$\leq 0.1\text{mg/kg}$
Phosphates	$\leq 0.5 \text{ mg/ I}$	$\leq 0.1\text{mg/kg}$
Conductivity	$\leq 15\mu\text{s /cm}$	$\leq 3 \mu\text{s /cm}$
PH Value	5 – 7.5	5 - 7
Appearance	Colorless, clean	Colorless, clean
Hardness	0.02 mmol/ I	0.02 mmol/ I





Appendix 2**Diagrams of the sterilization programs**

Programs Clave23+	Temperature (°C)	Pressure (kPa)	Holding time (min)	Total time (min)	Type	Max load (kg)	Max load per tray (kg)
SOLID (Unwrapped)	134	210	4	15-35	Unwrapped solid material	5.00	1.50
	121	110	20	31-40			
WRAPPED	134	210	6	30-50	Single-wrapped solid or hollow material	4.00	1.20
	121	110	20	48-60	Single-wrapped solid or hollow material	4.00	1.20
TEXTILE	134	210	8	36-55	Unwrapped porous material	1.25	0.40
					Single-wrapped porous material	1.10	0.30
	121	110	30	55-70	Dual-wrapped porous material	0.75	0.25
					Single-wrapped hollow material	4.00	1.25
					Dual-wrapped solid and hollow material	2.00	0.60
PRION	134	210	18	50-75	Unwrapped porous material	1.25	0.40
					Single-wrapped porous material	1.10	0.30
	121	110	30	55-70	Dual-wrapped porous material	0.75	0.25
					Single-wrapped hollow material	4.00	1.25
					Dual-wrapped solid and hollow material	2.00	0.60
LIQUID (optional)	134	210	10	35-60	Liquid	1.20	0.40
	121	110	30	35-65			
Drying (optional)	—	—	—	1-20	—	—	—

The time required for sterilizer to be ready for routine use after the power is switched on is less than 15 minutes.


The maximum temperature of the 134°C sterilization cycle is 137°C

The maximum temperature of the 121°C sterilization cycle is 124°C



Appendix 3

Wiring Diagram




TP1: Steam generator temperature sensor
 TP2: Inner temperature sensor
 TP3: Temperature sensor of chamber wall

V1: Vacuum pump valve
 V4: Water release valve

Appendix 4

Hydraulic Diagram



V1: Air release valve (normally open) P/N 90012077
 V4: Water release valve (normally closed) P/N 90012025



Alternate V1, p/n 90012024



Air release valve
 V1



Water release valve
 V4



Certificate of Warranty

All new Clave16 or Clave 23 or Class B autoclaves installed by a Flight Dental Systems authorized dealer for a period of two (2) full year from the time of purchase. The warranty covers defects in parts, workmanship and materials for two (2) year except door gaskets and filters which are wear and tear items. This warranty does not include labor or installation. This warranty does not apply to any device that has been subjected to improper use or accident; nor shall it extend to autoclaves that have been repaired or altered by an unauthorized Flight dealer or technician. The warranty also does not include routine cleaning or preventive maintenance.

Flight's obligation is limited to the repair or replacement of parts for the autoclave. No other warranties or obligations are expressed or implied. The user must follow the instructions for use as outlined in the user manual. To activate the warranty, the registration card must be completed and mailed or faxed to Flight within fourteen (14) days of purchase or you may call our customer service department at the number listed below. Products will only be received and accepted for repair from an authorized dealer and only with prior return authorization from Flight.

All Transportation charges to and from Flight must be paid by the owner of the Autoclave. Flight will not accept COD shipments. If repairs are needed during the first 90 days after purchase of their autoclave and a local authorized service dealer is not available. Flight will arrange pick up of the unit at Flight's expense. This will be on an individually evaluated basis and ONLY with pre-approval. Note: If you have any questions or there are any difficulties with this instrument and the solution is not covered in this manual, please contact your dealer or Flight Dental Systems. Do no attempt to service this device yourself.

Authorized Dealer: _____

Installed by: _____

Product Serial Number: _____

Product Description: _____

Product Model: _____

Purchased Date: _____ Invoice Number: _____

End User Name: _____

Telephone: _____ Fax: _____

Email: _____

Address: _____

City: _____ State/Province: _____

Zipcode: _____ Country: _____

A division of HR Dental Products
21 Kenview Blvd, Unit 11, Brampton, Ontario, Canada, L6T 5G7
Tel: 905-799-0517- Toll Free N.A.: 1-866-799-0517- Fax: 905-793-2406

Ver.1.2- 11/10/17

Contact Information

Thank you for taking the time to review the Flight Dental Unit's Instruction Manual. Your feedback or comments regarding the document is welcomed. For any comments and concerns please mail, e-mail or phone us at:

Sales and Marketing

Joseph Hui
1-866-799-0517
905-799-0517

Technical Support

1-866-799-0517
905-799-0517

Parts and Warranty

1-866-799-0517
905-799-0517

Mailing and Shipping Address (Office)

21 Kenview Dr, Unit 11
Brampton, ON Canada L6T 5G7
L6T 5G7