
	CAREGLOVE GLOBAL SDN BHD (933760-W)	No. / Rev.	SC-598/00
	PRODUCT SPECIFICATION		Page

240MM, 3.5G NITRILE EXAMINATION, POWDER FREE, FINGERTIP TEXTURED GLOVE, ND2DPFT1B2

DESIGN & FEATURE					
Product Code	ND2DPFT1				
Material Type	Nitrile				
Color	Blue				
Glove Type	Powder Free				
Former Type	Fingertip Textured				
Hand Design	Ambidextrous				
PHYSICAL DIMENSION					
Parameter	Size	EN-455-2 (2015)		TYPICAL VALUE	
Palm Thickness, mm (Single layer)	All sizes	-		0.05 – 0.09	
Finger Thickness, mm (Single layer)	All sizes	-		0.07 – 0.11	
Palm Width, mm	XS	70 ± 10		70 ± 5	
	S	80 ± 10		85 ± 5	
	M	95 ± 10		95 ± 5	
	L	110 ± 10		105 ± 5	
Length, mm	XS	Min 240		Min. 240	
	S	Min 240		Min. 240	
	M	Min 240		Min. 240	
	L	Min 240		Min. 240	
PHYSICAL PROPERTIES					
		EN-455-2 (2015)		TYPICAL VALUE	
Before Aging					
Force at Break, N	All sizes	Min. 6		Median. 6	
After Aging					
Force at Break, N	All sizes	Min. 6		Median. 6	
RESIDUAL POWDER CONTENT (ASTM D6124)					
Residual Powder Content, mg/glove	All sizes	Max. 2		Max. 2	
PRESHIPMENT INSPECTION		EN-455-1 (2000)		CGG SPECIFICATION	
		Inspection Level	AQL	Inspection Level	AQL
Barrier Defect (Water Leak Test)	All sizes	G-1	1.5	G-1	1.5
Visual Defect - Major	All sizes	N/A	N/A	G-1	2.5
Visual Defect - Minor	All sizes	N/A	N/A	G-1	4.0
Physical Dimensions	All sizes	N/A	N/A	S-2	4.0
Physical Properties	All sizes	N/A	N/A	S-2	4.0
Powder Free Residue	All sizes	N/A	N/A	N=5	N/A

	CAREGLOVE GLOBAL SDN BHD (933760-W)	No. / Rev.	SC-598/00
	PRODUCT SPECIFICATION	Page	2 of 2

240MM, 3.5G NITRILE EXAMINATION, POWDER FREE, FINGERTIP TEXTURE GLOVE, ND2DPFT1B2

PRODUCT LABELING INFORMATION (IN-HOUSE)		
Lot Numbering	All sizes	In-house or per customer requirement
Expiry Date (if applicable)	All sizes	In-house or per customer requirement
Manufacture Date	All sizes	In-house or per customer requirement

Prepared by,
 Name : Farhana Aqilah Zulkipli
 Position : Quality System Executive
 Date : 10th October 2018

Reviewed by,
 Name : Siti Khatijah Ahmad
 Position: Quality Assurance Manager
 Date :10th October 2018