

Technical Data Sheet
Procedural Face Mask Meet ASTM F2100-19 Level 3 Requirement


Intended Use : Procedure masks are intended protect and reduce the risk of the transfer of microorganisms, body fluids and particulate material transfer. Latex free.

| Characteristic | Testing Method | Limit |
|------------------------------|-----------------------|--------------|
| Bacterial Filtration (BFE) | ASTM F2101 | ≥ 98% |
| Particulate Filtration (PFE) | ASTM F2299 | ≥ 98% |
| Fluid Resistance | ASTM F1862 | 160 mmHg |
| Differential Pressure | EN 14683 | < 6.0 |
| Flammability | 16 CFR 1610 | Class I |

Material:

| Mask Construction | Material |
|----------------------------|---|
| Outer Layer: | Fluid repelling non-woven spun bond polypropylene |
| Middle Layer (Filtration): | Non-woven melt-blown media |
| Inner Layer: | Non-woven spun bond polypropylene |
| Nose clip: | Malleable wire |
| Ear-loop: | Polyester |

Mask Dimension and visual:



Packaging information:



Picture for reference only

Mask Color: White

Box Quantity: 50 pcs per box

Carton Quantity: 20 boxes per carton

Special Instruction: Store under cool, clean, and dry conditions. Avoid excessive heat.


| Box Dimension | Measurements |
|----------------------|---------------------|
| Length | 210 mm / 8.26" |
| Width | 120 mm / 4.72" |
| Height | 90 mm / 3.54" |

| Carton Dimension | Measurements |
|-------------------------|---------------------|
| Length | 445mm / 17.5" |
| Width | 254mm / 10" |
| Height | 470mm / 18.5" |


Instructions:

Disposable mask is for single usage only


FRONT - face out



Metal nose band at the **TOP**
Pleats fold **DOWN** on the **OUTSIDE**



BACK - face in



Metal nose band at the **TOP**
Pleats fold **UP** on the **INSIDE**