



AdvanTech[®]

Sterilizers



1ST
EDITION



Thank you for your interest in AdvanTech Sterilizers. We are a leader in sterilization products and services, where our trademark name represents quality and excellence.

Our products and services make us valuable partners to health care providers, as we continue enhancing the value we offer in every aspect of our company.

Customers may contact our representatives for technical or sales assistance from Monday-Friday 6am-3:30pm PST.

We guarantee your 100% satisfaction on all our products and services. Join the thousands that rely on AdvanTech Sterilizers and contact us today!

CONTACT US

CALL TOLL-FREE: 800.701.0186
EMAIL: info@advantechsterilizers.ca
MAIL: 177-33771 George Ferguson Way
Abbotsford, BC V2S 6H1

BUY ONLINE www.AdvanTechSterilizers.ca

Your Single-Source Packaging Solution!



NOTE: Prices are subject to change at any time. See our website for our most up-to-date pricing!

**Crepe Paper &
SMS Wrap**
pages 4-6

**Indicator
Tape**
page 7

**Pouches/Bags
& Tubing**
pages 8-9

**Label Dispensers
& Refills**
pages 10-11

DID YOU KNOW?

14 billion

pounds (68 kg) of waste are dumped into the ocean every year. Most of it is plastic.

200 million

tons of plastic are produced by manufacturers per year. Less than 3.5 is recycled, meaning every year we add 193 million tons of plastic to the world - permanently.

35,500 days

for a SMS non-woven cassette wrap to fully decompose. How many pieces of SMS wrap do you use each day?

**How can YOU
make a difference?**

make a start by choosing the right wrap >>


Why Crepe Paper?

- Biodegradable product – FULLY decomposes in only 10 weeks while SMS Wrap takes over 100 years to decompose!
- Easy to detect compromised packaging! Holes or punctures show much clearer on crepe, compared to SMS wrap.
- More cost-effective method of wrapping instruments.
- Minimizes the carbon impact on the environment (carbon compensated for CO2 emissions)
- Renewable and sustainable (Crepe paper is made of cellulose which comes from 30-year old trees. Suppliers of cellulose replace one tree for every tree cut down.)
- Compatible with Steam, EO Gas and Gamma Radiation sterilization processes


Easily Detect Compromised Packs:

A wrapped tray/pack requires uncompromised packaging. A pack wrapped with crepe paper that has been punctured can be easily detected compared to SMS wrap (see pictures below). SMS wrap has a hole-pattern look, making it difficult to detect a pin hole.

Hole punctures in SMS wrap are unnoticeable




A close-up view of the hole punch in SMS wrap




Rapidly Biodegrades:

Once crepe paper is discarded, it easily biodegrades. Within 10 weeks, a sheet of crepe wrap is completely disintegrated! On the other hand, SMS wrap takes over 100 years to decompose!



CREPE AT START



AFTER ONE WEEK
OF COMPOSTING



AFTER THREE WEEKS
OF COMPOSTING



AFTER SIX WEEKS
OF COMPOSTING



AFTER TEN WEEKS
OF COMPOSTING



AdvanCheckTM Crepe Paper

AdvanCheck TM Crepe Paper			1-4 BOXES	5+ BOXES
W-CRP12	12" x 12" Crepe Wrap	1,000/bx	\$94.00/bx	\$89.00/bx
W-CRP15	15" x 15" Crepe Wrap	1,000/bx	\$119.00/bx	\$114.00/bx
W-CRP18	18" x 18" Crepe Wrap	500/bx	\$99.00/bx	\$94.00/bx
W-CRP20	20" x 20" Crepe Wrap	500/bx	\$101.00/bx	\$96.00/bx
W-CRP24	24" x 24" Crepe Wrap	500/bx	\$129.00/bx	\$124.00/bx

Prefer a crisp, thin paper to wrap your cassettes? Our standard crepe paper is the right choice for you. Request your free samples today! Please read the preceding page to see how you can make a difference.



AdvanCheckTM PREMIUM Crepe Paper

AdvanCheck TM PREMIUM Crepe Paper			1-4 BOXES	5+ BOXES
W-CRP12-P	12" x 12" PREMIUM Crepe	1,000/bx	\$114.00/bx	\$109.00/bx
W-CRP15-P	15" x 15" PREMIUM Crepe	1,000/bx	\$179.00/bx	\$174.00/bx
W-CRP18-P	18" x 18" PREMIUM Crepe	500/bx	\$134.00/bx	\$129.00/bx
W-CRP20-P	20" x 20" PREMIUM Crepe	500/bx	\$149.00/bx	\$144.00/bx
W-CRP24-P	24" x 24" PREMIUM Crepe	500/bx	\$199.00/bx	\$194.00/bx

Try our thicker, softer crepe paper that almost feels like cloth! PREMIUM Crepe provides a strong tear-resistant barrier of protection.

Please read the preceding page to see how you can make a difference.



AdvanCheckTM ComboWrap

AdvanCheck TM ComboWrap			1-4 BOXES	5+ BOXES
W-CW12	12"x12" ComboWrap	400 interleaved sheets/bx (200 sets)	\$115.00/bx	\$110.00/bx
W-CW15	15"x15" ComboWrap	400 interleaved sheets/bx (200 sets)	\$118.00/bx	\$113.00/bx
W-CW18	18"x18" ComboWrap	400 interleaved sheets/bx (200 sets)	\$122.00/bx	\$117.00/bx
W-CW20	20"x20" ComboWrap	300 interleaved sheets/bx (150 sets)	\$106.00/bx	\$101.00/bx
W-CW24	24"x24" ComboWrap	400 interleaved sheets/bx (200 sets)	\$175.00/bx	\$170.00/bx


The ComboWrap Solution:

INNERWRAP: Environmentally-friendly crepe paper for the first layer of wrap which also provides a superior barrier for greater protection of your instruments.

OUTERWRAP: A strong, fluid-resistant, breathable & durable protection of SMS for the second layer of wrap.

Why ComboWrap?

- Minimize your impact on our environment! Crepe paper biodegrades in 10 weeks (vs. 100 years for SMS wrap).
- SMS provides the strong, tear-resistant outer layer
- SMS and Crepe paper are interleaved throughout the box.



SMS Wrap
(outer layer)


Crepe Premium Paper
(inner layer)

CLICK HERE
To Request Wrap Samples!

Need help choosing the right wrap for you?
Our staff is here to help! Give us a call at 800.701.0186

AdvanCheck SMS Autoclave Wrap



AdvanCheck  Blue SMS Wrap			1-4 BOXES	5+ BOXES
W-SMS12	SMS Wrap, 12" x 12"	500 pc/bx	\$89.00/bx	\$84.00/bx
W-SMS15	SMS Wrap, 15" x 15"	500 pc/bx	\$104.00/bx	\$99.00/bx
W-SMS18	SMS Wrap, 18" x 18"	200 pc/bx	\$84.00/bx	\$79.00/bx
W-SMS20	SMS Wrap, 20" x 20"	200 pc/bx	\$86.00/bx	\$81.00/bx
W-SMS24	SMS Wrap, 24" x 24"	250 pc/bx	\$99.00/bx	\$94.00/bx
W-SMS30	SMS Wrap, 30" x 30"	200 pc/bx	\$114.00/bx	\$109.00/bx
W-SMS36	SMS Wrap, 36" x 36"	100 pc/bx	\$94.00/bx	\$89.00/bx

KEY FEATURES:

- High quality and consistent
- Strong, fluid-resistant, breathable and provides durable protection
- Manufactured by a world leading production facility
- Compatible with Steam, EO Gas and Gamma Radiation sterilization processes

CLICK HERE

To Request Wrap Samples!

KimGuard One-Step* Quick Check* Sterilization Wrap

OUR 2-IN-1 WRAP!




*one application =
two sheets of
KIM-GUARD* Sterilization Wrap

KimGuard One-Step* Quick Check* Wrap			1-4 CASES	5+ CASES
W-KC100-12	12" x 12" One-Step Wrap	480 sets/case	\$99.00/cs	\$94.00/cs
W-KC100-15	15" x 15" One-Step Wrap	480 sets/case	\$129.00/cs	\$124.00/cs
W-KC100-18	18" x 18" One-Step Wrap	480 sets/case	\$209.00/cs	\$204.00/cs
W-KC100-20	20" x 20" One-Step Wrap	480 sets/case	\$219.00/cs	\$214.00/cs
W-KC100-24	24" x 24" One-Step Wrap	240 sets/case	\$159.00/cs	\$154.00/cs
W-KC100-30	30" x 30" One-Step Wrap	144 sets/case	\$189.00/cs	\$184.00/cs

KEY FEATURES:

- Confidence you can see, quicker than ever. Quick Check wrap is the new fast, easy way to ensure instrument sterility is not compromised. Two colours – white and blue – bonded together, makes inspecting trays is easier and faster than ever!
- Wrap was invented and is still fully manufactured in the USA
- Resistance to Linting
- Wraps Up to 3 lbs / 1.36 kgs (other grades available)

Sterilizer Indicator Tape Dispenser



Sterilizer Indicator Tape Dispenser		1 EACH
T-DISP2	Tape Dispenser, stores up to 2 rolls of 1" tape per dispenser	\$65.00/ea
T-DISP3	Tape Dispenser, stores up to 3 rolls of 1" tape per dispenser	\$69.00/ea
T-DISP4	Tape Dispenser, stores up to 4 rolls of 1" tape per dispenser	\$79.00/ea

Tape dispenser shown at right can hold between two and four rolls of indicator tape at one time to suit your needs. Each roller holds tape as wide as 1".



Dispenser Can Hold
2-4 Rolls of Tape!



Blade Protector
If Not In Use

SHARP BLADE

Steam Indicator Tape

			ORDER BY ROLL	CASE PRICE!	
Sterilizer Indicator Tape			1-9 ROLLS	10+ ROLLS	1+ CASE
T-S.5	½" x 180' Tape for Steam	72 rolls/cs	\$6.50/roll	\$6.00/roll	\$5.50/roll
T-S.75	¾" x 180' Tape for Steam	48 rolls/cs	\$8.50/roll	\$8.00/roll	\$7.50/roll
T-S1	1" x 180' Tape for Steam	36 rolls/cs	\$9.50/roll	\$9.00/roll	\$8.50/roll

Indicator tapes for Steam with superior adhesion to woven and non-woven fabrics. Tape adheres to items throughout the sterilization process. Steam tape indicator turns from white-to-brown/black when exposed to steam sterilization conditions.



STOCK UP & SAVE
take advantage of
quantity pricing!

Green Lead-Free Indicator Tape

			ORDER BY ROLL	CASE PRICE!	
Green Lead-Free Indicator Tape			1-9 ROLLS	10+ ROLLS	1+ CASE
T-S.75LF	Lead Free ¾" x 180' Tape for Steam	48 rolls/cs	\$9.50/roll	\$9.00/roll	\$8.50/roll
T-S1LF	Lead Free 1" x 180' Tape for Steam	36 rolls/cs	\$10.50/roll	\$10.00/roll	\$9.50/roll



Our environmentally-friendly green tape does not contain lead. Tape adheres to items throughout the sterilization process. Designed for gravity sterilizers at 121°C for 30 minutes and pre-vacuum steam sterilizers at 132°C for 3 minutes and 135°C for 3 minutes.



Professionals Choice Steam Sterilizer Pouches/Bags



Exterior & Interior
Class 4 Chemical Indicators!

CLICK HERE
To Request a Sample!

			BOX PRICE!	CASE PRICE!
Class 4 Steam Pouches			PRICE/BOX	1+ CASES
P4-442089	2¼" x 4"	200/bx, 12 boxes/cs	\$9.50/bx	\$8.50/bx
P4-442090	2¾" x 10"	200/bx, 6 boxes/cs	\$11.50/bx	\$10.50/bx
P4-442100	3½" x 6½"	200/bx, 6 boxes/cs	\$11.75/bx	\$10.75/bx
P4-442101	3½" x 10"	200/bx, 6 boxes/cs	\$12.50/bx	\$11.50/bx
P4-442102	5¼" x 11"	200/bx, 6 boxes/cs	\$15.25/bx	\$14.25/bx
P4-442103	7½" x 13"	200/bx, 5 boxes/cs	\$29.00/bx	\$28.00/bx

KEY FEATURES:

- Meets infection control guidelines with interior and exterior Class 4 indicators for steam sterilization
- Manufactured in Canada
- Lead-free process indicators
- Superior medical-grade paper repels tearing
- Easy peel-and-seal with perforated folding tab


Class 4 Internal/External multi-parameter indicators change colour once all three parameters for sterilization have been met:
TIME, TEMPERATURE & PRESSURE

MADE IN CANADA

Class 1 Steam Sterilizer Pouches/Bags

			BOX PRICE!	CASE PRICE!
Class 1 Steam Pouches			PRICE/BOX	1+ CASES
P1-SP0507	5¼" x 7½"	200/box, 12 boxes/cs	\$15.00/box	\$12.00/bx
P1-SP1015	10" x 15"	200/box, 8 boxes/cs	\$62.00/box	\$49.50/bx
P1-SP1215	12" x 15"	200/box, 5 boxes/cs	\$72.00/box	\$65.00/bx
P1-SP1218	12" x 18"	100/box, 10 boxes/cs	\$65.00/box	\$57.00/bx

Sterilization Tubing



			ROLL PRICE!	CASE PRICE!
Sterilization Tubing			PRICE/ROLL	1+ CASES
ST-S2	2" x 100' Tubing, 10 rolls/cs		\$11.95/roll	\$9.50/roll
ST-S3	3" x 100' Tubing, 10 rolls/cs		\$12.95/roll	\$10.50/roll
ST-S4	4" x 100' Tubing, 10 rolls/cs		\$13.95/roll	\$11.50/roll
ST-S6	6" x 100' Tubing, 10 rolls/cs		\$17.95/roll	\$15.50/roll
ST-S9	9" x 100' Tubing, 10 rolls/cs		\$20.95/roll	\$18.50/roll

Tubing is transparent on the top side allowing you to identify items in each package, and superior medical-grade paper that allows for sterilization penetration. Our high quality tubing comes with EO Gas and Steam Indicators for the sterilization process.

Sterilization Tubing Heat Sealer



Sterilization Tubing Heat Sealer			1 EACH
ST-HS	Heat Sealer, 110V, 8" (20cm) wide, 1 each		\$349.00/ea

This sealer creates an air and watertight seal on sterilization rolls (paper-plastic) and sterilization tubing (plastic-plastic) materials. This compact unit has a built-in warning alarm which notified the user when each seal is complete. The seal produced is an embossed, honeycomb texture ⅜" (10mm) wide. This medical impulse sealer is designed especially for the Health Care Professionals, that require wider and stronger seals.



Dry Heat Sterilizer Pouches

Class 1 Dry Heat Sterilization Pouches		1-4 BOXES	5+ BOXES
P1-DH210	2" x 10" Nylon Pouch, 100/box	\$57.00/bx	\$54.00/bx
P1-DH310	3" x 10" Nylon Pouch, 100/box	\$62.00/bx	\$59.00/bx
P1-DH410	4" x 10" Nylon Pouch, 100/box	\$66.00/bx	\$63.00/bx
P1-DH710	7" x 10½" Nylon Pouch, 100/box	\$82.00/bx	\$79.00/bx
P1-DH913	9½" x 13" Nylon Pouch, 100/box	\$99.00/bx	\$96.00/bx

The self-seal clear nylon pouches each include Class 1 chemical indicators. Easy peel-and-seal folding tab.



Steri-SURE™ Dry Heat Sterilizer Bags

Class 1 Dry Heat Bags		1-4 BOXES	5+ BOXES
P1-DHB2	2½" x 1½" x 10½" Sterilizer Gusset Bg, 250/bx	\$83.00/bx	\$80.00/bx

The Steri-SURE™ paper sterilization bags each include Class 1 chemical indicators.



Dry Heat Sterilization Tubing

Dry Heat Nylon Sterilization Tubing		1-4 ROLLS	5+ ROLLS
ST-DH2	2"x100' (5cm x 30M) Tubing	\$10.99/roll	\$10.49/roll
ST-DH3	3"x100' (7.5cm x 30M) Tubing	\$13.99/roll	\$13.49/roll
ST-DH4	4"x100' (10cm x 30M) Tubing	\$17.99/roll	\$17.49/roll
ST-DH6	6"x100' (15cm x 30M) Tubing	\$25.99/roll	\$25.49/roll

A clear, strong, puncture resistant, reusable autoclave tubing available in four different widths that can be cut to accept any instrument length. A strong, reliable seal is formed when used in conjunction with a heat sealer (see page 8). Nylon tubing is available specifically for use in Dry Heat and Chemical Vapour sterilizers.



For Dry Heat &
Chemical Vapour
Sterilization

Sterilization Tubing Dispenser

Sterilization Tubing Dispenser		1 EACH
ST-DR-6	Tubing Dispenser, 1 ea (Holds 2 rolls, and up to 6" wide)	\$59.00/ea
ST-DR-10	Tubing Dispenser, 1 ea (Holds 3 rolls, and up to 10.5" wide)	\$69.00/ea

This dispenser fits all sterilization tubing core sizes, from 1" - 3" in diameter. The suction-cup feet secure the dispenser to the table.



Labelex™ Dispensers

ASK ABOUT OUR DISPENSER TRADE-IN PROGRAM!

Receive **50% OFF** any Labelex dispenser by trading in your used dispenser of another brand! Ask us how today!

SAVE ON YOUR LABEL REFILLS!
AS LOW AS \$9.95/ROLL

Label Dispenser		1-2 DISPENSERS	3+ DISPENSERS
L-X20-R	RED Labelex Dispenser	\$235.00/ea	\$215.00/ea
L-X20-B	BLUE Labelex Dispenser	\$235.00/ea	\$215.00/ea
L-X20-G	GREEN Labelex Dispenser	\$235.00/ea	\$215.00/ea
L-X20-IR	Replacement Ink Rolls for Labelex Disp., 5/pk	\$28.00/pk	-
L-Z20-IR	Replacement Ink Rolls for Black Dispenser, 5/pk	\$28.00/pk	

The perfect system to organize and track your sterilized bags. Great for any size office. Try out our user-friendly, long-lasting, high quality label dispenser today!



Labelex™ Benefits & Features



EFFORTLESS

Save the tediousness of handwriting the date, time, load number and tray number on every pouch. With a click of the button your job is done!



SAVE YOUR PRACTICE TIME

No more handwriting your information onto your packs! Save yourself and your staff time by automating your practice!



TRACEABILITY

If a sterilization error occurs, sterilized packs or pouches can be easily traced, from the sterilizer to the patient's records!




EASY TO LOAD

This dispenser is notably easier to load than other dispenser brands. We also include a user-friendly instruction guide with your order.



QUALITY

 Labelex™ dispensers are built to last. All dispensers and labels are made in Italy. Strong label adhesive keeps label in place during and after sterilization.



GUARANTEED

All label dispensers come with a lifetime warranty* and will be replaced at no cost to you. We stand behind our quality products!



LOW COST

Our label dispensers are reasonably priced and are considerably lower than most leading dispenser brands. The dispenser refills are cost-efficient; as low as **\$9.95 for 1,000 labels!**



EASY IDENTIFICATION

Easy to identify which dispenser belongs to each sterilizer!

NEW

Tip from the pros: If your office has more than one sterilizer, use a dispenser for each sterilizer, with coordinating coloured labels.

Single-Ply Labelex Refills

Single-Ply Labels		1-3 PACKS (\$10.95/ROLL)	4-6 PACKS (\$10.45/ROLL)	7+ PACKS (\$9.95/ROLL)
L-Z20-SA	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack
L-Z20-SU	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack
L-Z20-SS-RD	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack
L-Z20-SS-BL	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack
L-Z20-SS-GR	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack
L-Z20-SS-WH	1,000/roll, 12 rolls/pack	\$131.40/pack	\$125.40/pack	\$119.40/pack

Designed for computerized record keeping!

Apply the label to the package before or after sterilizing to track the date of sterilization and the load number.

L-Z20-SU

L-Z20-SA

L-Z20-SS-BL
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged, Opened or Wet

L-Z20-SS-GR
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged, Opened or Wet

L-Z20-SS-RD
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged or Opened

L-Z20-SS-WH

~ Other options available ~
Speak to one of our representatives today!

Dual-Ply Labelex Refills

Dual-Ply Labels		1-3 PACKS (\$11.45/ROLL)	4-6 PACKS (\$10.95/ROLL)	7+ PACKS (\$10.45/ROLL)
L-Z20-DA	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack
L-Z20-DU	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack
L-Z20-DS-RD	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack
L-Z20-DS-BL	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack
L-Z20-DS-GR	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack
L-Z20-DS-WH	500/roll, 12 rolls/pack	\$137.40/pack	\$131.40/pack	\$125.40/pack

Designed for paper record keeping!

Apply the label to the package to be sterilized. At time of use, remove the label from the wrap or bag and place the label onto the patient's record.

L-Z20-DU

L-Z20-DA

L-Z20-DS-BL
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged, Opened or Wet


L-Z20-DS-GR
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged, Opened or Wet

L-Z20-DS-RD
DATE / LOAD NO.
STERILIZED
Indefinite Shelf Life
Unless Damaged or Opened

L-Z20-DS-WH


~ Other options available ~
Speak to one of our representatives today!

- Dual-Ply Label Advantage -



Heli-Check Hollow Load PCD

As low as **\$.79/test!**



Heli-Check Hollow Load PCD		1 KIT	2 KITS	3+ KITS
PCD-HC	Helix Test Device with 250 Heli-Check tests	\$259.00/kit	\$249.00/kit	\$239.00/kit

A hollow load Process Challenge Device for table-top pre-vacuum steam sterilizers. A reusable device that can be utilized in each sterilization load as an independent control device for decisions regarding load release.

- Verifies steam penetration and exposure levels in lumens
- Re-usable test device - can be used up to 250 times
- Low cost per test
- Results can be stored as a permanent record.
- Helix devices are not useful for testing of porous loads and therefore should not be used as a Bowie-Dick test
- Conforms to EN 867-5 and ISO 11140-1, Class 2

INTERPRETATION



The presence of any pink colour remaining indicates insufficient sterilization conditions and re-sterilization is recommended.

- Load release common failures:
- Inadequate cycle time
 - Improper cycle selection
 - Overloading the sterilizer
 - Poor steam quality
 - Equipment malfunction

What is a Process Challenge Device (PCD)?

Process challenge device (PCD): An item designed to simulate product to be sterilized and to constitute a defined challenge to the sterilization process and used to assess the effective performance of the process. A PCD is a challenge test pack or test tray that contains a biological indicator, a Class 5 integrating indicator, or an enzyme-only indicator.

Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008



Are you monitoring your cleaning process?



Detergents &
Ultrasonic Cleaners
pages 15-17

HLD Sterilants/Monitoring
& Instrument Care
pages 18-19

Washer Disinfector
Monitoring
pages 20-21

INSTRUMENT CARE USING EMPOWER AND METRILUBE PRECLEAN, LUBRICATE AND MONITOR YOUR HIGH-LEVEL DISINFECT/STERILIZE

Preclean...

using EmPower® Foam™

1. Place items in clearly marked decontamination area.
2. Thoroughly spray EmPower® Foam™ directly onto instruments to be cleaned.
3. Allow instruments to remain covered with enzymatic foam to dissolve bioburden.
4. Reapply as needed to keep instruments moist.
5. Rinse instruments thoroughly.

using EmPower® Enzymatic Detergent

1. Place items in clearly marked decontamination area.
2. Dilute EmPower® Enzymatic Detergent 1 fluid ounce (1 pump) per gallon of water.
3. Immerse instruments until soil is dissolved and removed. If matter is dried on, soak longer until visibly clean.
4. Manually clean with brush or irrigate lumens under the solution level.
5. Rinse instruments thoroughly.
6. Air-dry or dry with clean towel.
7. Inspect instruments. If instruments are not clean, repeat steps.



Spray EmPower® Foam™ to conveniently remove organic debris while reducing splashes and messy spills.

Lubricate...

using MetriLube® Instrument Lubricant Ready-to-Use/Concentrate

1. Apply MetriLube® directly onto clean instruments.
2. Wipe dry with clean towel.



Use EmPower® enzymatic detergent to effectively remove organic debris in both manual and automated reprocessing systems.

High-Level Disinfect/Sterilize...

using MetriCide® and MetriTest™ High-Level Disinfectant/Sterilant

1. Prepare instruments for sterilization or high-level disinfection.
2. Sterilize or high-level disinfect according to the manufacturer's labels.

Verify/Monitor...

using DISINTEK™ and MetriTest™ monitoring strips

1. Using monitoring strips specific to each solution. Dip the strip into the soaking tray to ensure that the concentration is at the minimum effective concentration (MEC) for that solution.
2. Follow manufacturer's instructions to see how long to hold strip in solution.
3. 60-90 seconds later, compare colour of indicator pad to the colour chart.



Use MetriLube® Spray or concentrate for lubricating instruments to help prevent spotting, staining, rusting and freezing up.

Storage...

Store instruments according to your healthcare system's policies and regulations.

Note: Follow standard precautions for personal protection as appropriate.

EmPower® Foam™ Enzymatic Spray

		ORDER BY BOTTLE			ORDER BY CASE
EmPower® Foam™ Spray		1-3 BOTTLES	4-6 BOTTLES	7-11 BOTTLES	1 CASE
CS-EPF-24	EmPower® Foam™ - 24 oz. Spray, 1 bottle, 12/cs, (MFR# 10-4224)	\$33.00/btl	\$32.00/btl	\$31.00/btl	\$29.00/btl

EmPower® Foam™ propels thick foam of protease enzymes. Effectively cleans away blood, tissue, mucous and other protein-rich body fluids. Foaming spray is ideal for use in surgery, central processing, operating room and other departments where instrument cleaning is not immediately available.

KEY FEATURES & BENEFITS:

- Ideal for instruments with presence of moderate to heavy organic debris
- Designed for manual reprocessing
- Ready to use in a convenient sprayer
- Thick foam reduces splashing and spills
- Covers instruments to quickly break down bio-burden



EmPower® Dual-Enzymatic Detergent

		ORDER BY BOTTLE	ORDER BY THE CASE & SAVE!		
EmPower® Detergent		1 GALLON	1 CASE	2 CASES	3+ CASES
CS-EPS-G	EmPower® Solution, 1 gallon, 4 gallons/case (MFR# 10-4100)	\$56.00/gal	\$55.00/gal	\$54.00/gal	\$53.00/gal

EmPower® is a fresh scented, detergent cleaner containing two protease enzymes that offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments. Non-corrosive to metals.

KEY FEATURES & BENEFITS:

- Dual-enzymatic detergent cleaner
- Low-foaming formula is ideal for use in automated machines
- Less corrosive than distilled water
- Each case includes one pump nozzle that dispenses 1 oz. of fluid each pump.
- One gallon produces 128 gallons of solution (1 oz. of empower, to each gallon of distilled water)



Ultrasonic Cleaner Tablets

Tuttnauer Ultrasonic Tablets		1 BOX	2 BOXES	3+ BOXES
UC-CST	Clean&Simple Ultrasonic Cleaner Solution, 64 Tablets/bx	\$105.00/bx	\$100.00/bx	\$95.00/bx



Clean & Simple is an ultrasonic/enzymatic solution dispensed in a unique tablet form. Drop 2 tablets into your ultrasonic tank already filled with tap water and the tablets will dissolve in about a minute. This is a non-ionic solution which is non-corrosive, non-toxic, biodegradable and non-caustic. No bulky storage or pre-mixing required. Eliminates handling of heavy 1 gallon containers.

Economical Ultrasonic Cleaner



Special Price!
ONLY \$199.00

1.4 Litre Tank
UC1.4L

Removable Mesh Tray



Inside Unit



Ultrasonic Cleaners		Tank Size (LxWxH)	Watts	1 EACH
UC1.4L	1.4L Tank (1.5 Quarts)	6½" x 5½" x 2½"	60W	\$199.00/ea

Our 1.4L Digital Ultrasonic Cleaner system uses the process of "cavitation" in cleaning items which is more efficient and have much shorter cleaning time.

KEY FEATURES INCLUDE:

- Comparable with other brands
- 1 Year Warranty
- Built-in programmer that easily allows you to adjust 5 preset cycles
- Super quiet, solid-state circuitry
- Isolated-circuit design keeps water away from circuit for enhanced safety
- Removable mesh tray to prevent damage as well as positioning the item within the optimal cleaning zone of the tank.
- High power transducer for superior cleaning results
- 1.47 quart (1.4L) chamber

Premium Ultrasonic Cleaners



Ultrasonic Cleaners		Tank Size (LxWxH)	Watts	1 EACH
UC.7L	0.7L Tank (0.7 Quarts)	6" x 3½" x 2½"	50W	\$249.00/ea
UC2L	2L Tank (1.7 Quarts)	6" x 5¼" x 3¾"	380W	\$379.00/ea
UC3L	3L Tank (~.75 Gallons)	9" x 5½" x 4"	400W	\$579.00/ea
UC4L	4L Tank (~1 Gallon)	12" x 6" x 3¾"	750W	\$769.00/ea
UC6L	6L Tank (1.5 Gallons)	12" x 6" x 6"	780W	\$869.00/ea
UC10L	10L Tank (2.5 Gallons)	11" x 8" x 5"	600W	\$1,399.00/ea

KEY FEATURES INCLUDE:

- Free internal basket
- 2 year warranty
- Stamina XP Series
- Digital ultrasonic timer
- 110V, 40 Khz
- Improved quality, functionality and durability over other brands
- Engineered, designed, and manufactured in the U.S.A.
- Temperature regulated heating element for best cleaning results
- Reliable Industrial-Grade Transducer

Free Internal Basket



OTHER SIZES AVAILABLE!



.7 Litre Tank
UC.7L



2 Litre Tank
UC2L



3 Litre Tank
UC3L



4 Litre Tank
UC4L



6 Litre Tank
UC6L



10 Litre Tank
UC10L

WashChecks® U Monitors - *for your Ultrasonic Cleaner*

Wash-Checks® U Monitors		1 EACH	2 EACH	3+ EACH
WC-UC	50 Monitors/bag	\$85.00/bg	\$80.00/bg	\$75.00/bg
WC-H	Wash-Checks Holder, 1 Each	\$45.00/ea	\$42.00/ea	\$39.00/ea

Disposable ultrasonic cleaning monitors & reusable holders for monitoring the ultrasonic cleaning process.

Check the effectiveness of the wash process in your ultrasonic cleaner with a consistent monitor on a routine basis. When cleaning conditions are effective, test soil on monitor is dissolved.

Multi-parameter test

The presence of a large amount of blue test soil remaining on the monitor reflects a serious deficiency in the ultrasonic cleaning efficiency. It tests the effectiveness of time, temperature and the detergent.

Use with reusable Wash-Checks holder

Save money by using the same reusable holder for both the Wash-Checks® and Wash-Checks® U. Holder acts as a hinged instrument, measuring impingement.



Blue indicator will have completely dissolved once proper washing conditions have been met.

Are You Confident That Your Ultrasonic Cleaner Is Effectively Cleaning Your Instruments?

The WashChecks® U's are inserted in the reusable holder so that the blue square is fully exposed.

For Daily Machine Release: The monitor and holder are placed in the ultrasonic cleaner with cleaning liquid but no instruments with the blue square facing up, and a normal ultrasonic cycle is conducted.

For Load Release: The monitor and holder are placed in the cleaner on top of the instruments with the blue square facing up.



THE ONLY RELIABLE METHOD TO ENSURE THE ULTRASONIC CLEANER IS WORKING ADEQUATELY.

INTERPRETATION GUIDE

Wash-Checks U
Ultrasonic Monitor



UNUSED

Wash-Checks U
Ultrasonic Monitor



FAILED

Wash-Checks U
Ultrasonic Monitor



PASSED

MetriLube® Instrument Lubricant



		ORDER BY BOTTLE			ORDER BY CASE
MetriLube® Lubricant Spray		1-3 BOTTLES	4-6 BOTTLES	7-11 BOTTLES	1 CASE
CS-ML-24	MetriLube® Spray, 24 oz., 12 bottles/case (MFR# 10-3425)	\$24.00/btl	\$23.00/btl	\$22.00/btl	\$21.00/btl

MetriLube® is a concentrated oil-based instrument lubricant which coats metal instruments to prevent corrosion and help hinges work freely. Available in a 24 oz spray for manual use or in various gallon sizes for automated reprocessors.

KEY FEATURES & BENEFITS:

- Instrument lubricant concentrate
- Extends life of instruments
- Anti-corrosive
- Non-greasy and silicone-free
- Ideal for use in manual reprocessing or in automated systems
- One gallon produces 15 gallons of MetriLube spray

		ORDER BY BOTTLE	ORDER BY THE CASE & SAVE!		
MetriLube® Lubricant Refill		1 EACH	1 CASE	2 CASES	3+ CASES
CS-ML-G	MetriLube®, 1 gallon, 4 gallons/case (MFR# 10-3400)	\$56.00/gal	\$55.00/gal	\$54.00/gal	\$53.00/gal

Instrument Cleaning Brushes



		ORDER BY BRUSH	ORDER BY PACK
Nylon & Stainless Steel Brushes		1 BRUSH	PACK OF 3
CS-NB-Narrow	Nylon Cleaning Brush, 7" long	\$13.50/each	\$39.00/3 brushes
CS-SB-Narrow	Stainless Steel Brush, 7" long	\$12.00/each	\$35.00/3 brushes
CS-WB-Wide	Nylon Cleaning Brush, 3"x3" white	\$36.00/each	\$105.00/3 brushes

These brushes are intended to aid in the general cleaning of surgical instruments, devices and equipment. Brushes are latex free and require no special handling.

Nylon Bristles:

7" long, ideal for cleaning hard to reach places without scratching or harming delicate instrument surfaces. If desired, they may be sterilized.

Stainless Steel Bristles:

7" long, ideal for cleaning rasps, files, and burs. Bonus: Can be sterilized in a regular autoclave cycle!

White Nylon Bristles:

Surface cleaning brush: 3"x3" brush head, 4½" plastic handle, 2" white medical grade nylon bristles. For Use on sinks, autoclaves, ultrasonic cleaners, bowls, basins, and more!

High-Level 14-Day Disinfectant & Sterilant

		ORDER BY GALLON	ORDER BY THE CASE & SAVE!		
MetriCide® 14-Day Disinfectant		1 GALLON	1 CASE	2 CASES	3+ CASES
CS-MC14-G	14-Day MetriCide®, 1 gallon (3.8L), 4 gallons/case (MFR# 11-1401)	\$43.00/gal	\$42.00/gal	\$41.00/gal	\$40.00/gal

MetriCide® is a 2.6% buffered glutaraldehyde that can be used and reused as a sterilant and high-level disinfectant on immerseable, lensed and delicate instruments for up to 14 days.

KEY FEATURES & BENEFITS:

- Recommended for lensed instruments such as flexible and/or rigid scopes
- Compatible with heat-sensitive medical equipment and instrumentation
- Contains no surfactants
- No dilution required



Disintek™ GTA 1.5% Test Strips


Disintek™ 1.5% Indicator Tests		1 BOTTLE	2 BOTTLES	3+ BOTTLES
CS-1.5TS	1.5% Indicator Strips, 50/Bottle	\$116.00/btl	\$111.00/btl	\$105.00/btl

GTA 1.5% Test Strips are used to determine whether the glutaraldehyde concentration is at, above or below the Minimum Effective Concentration (MEC) of 1.5% for the following high-level disinfectants. 180 day shelf-life (double other brands). Cost-effective strips; use 70% less strips for quality control compared to other HLD tests!

FOR USE IN:

- CIDEX® Activated Dialdehyde Solution
- MetriCide® 14-day Sterilizing & Disinfecting Solution

PASS - the indicator pad is purple
FAIL - the indicator pad displays distinct white areas



Test Time: 90 Seconds



High-Level 28-Day Disinfectant & Sterilant

		ORDER BY GALLON	ORDER BY THE CASE & SAVE!		
MetriCide® 28-Day Disinfectant		1 GALLON	1 CASE	2 CASES	3+ CASES
CS-MC28-G	28-Day MetriCide®, 1 gallon (3.8L), 4 gallons/case (MFR# 11-2801)	\$33.00/gal	\$32.00/gal	\$31.00/gal	\$30.00/gal

MetriCide® 28 is a long-life, 2.5% buffered glutaraldehyde, that can be used and reused as a sterilant and high-level disinfectant for up to 28 days.

KEY FEATURES & BENEFITS:

- Recommended for immerseable and heat-sensitive equipment and instruments
- Contains a surfactant
- No dilution required



MetriTest™ GTA 1.8% Test Strips

MetriTest™ Indicator Tests		1 BOTTLE	2 BOTTLES	3+ BOTTLES
CS-1.8TS	MetriTest™ Indicator Strips, 60/Bottle (MFR# 10-304)	\$105.00/btl	\$100.00/btl	\$95.00/btl



MetriCide® 28 must maintain the minimum effective concentration (MEC) level of 1.8% glutaraldehyde or greater during the 28-day use life to effectively sterilize and high-level disinfect.

FOR USE IN:

- MetriCide® 28-Day Solution

PASS - the indicator pad is purple
FAIL - the indicator pad displays yellow areas



Test Time: 60 Seconds



DO YOU KNOW

If Your Washer-Disinfector
is Effectively Cleaning
Every Instrument?



WashChecks® Monitors - *for your Washer-Disinfector*

WashChecks® Monitors		1 EACH	2 EACH	3+ EACH
WC-WD	100 Monitors/bag	\$169.00/bg	\$159.00/bg	\$149.00/bg
WC-H	Wash-Checks Holder, 1 Each	\$45.00/ea	\$42.00/ea	\$39.00/ea

Disposable Wash Monitors & Reusable Holders for Monitoring the Wash Process in Washer-Disinfectors

Check your wash/detergent system effectiveness with a consistent monitor on a routine basis. If insufficient time, temperature, detergent concentration, spray arm function or enzyme soak occurs, a **RED** residue remains.

Test soil is designed to parallel the removal of blood and tissue from surgical instruments. No actual blood components from any source are contained in the test soil.

Use with reusable Wash-Checks® holder. Holder acts as a hinged instrument, measuring impingement.

Third party testing has proven Wash-Checks® to be the most reliable and effective wash monitors in the industry!



Red indicator will have completely dissolved once proper washing conditions have been met.

It LOOKS clean.
Is visual examination good enough?



Sometimes not...
And if it's not clean,
it's not sterile.



INTERPRETATION GUIDE

	UNUSED
	FAILED
	FAILED
	FAILED
	FAILED
	PASSED

The test soil (red dot) simulates human blood and tissue debris. When the test is placed in the holder the test soil is half obscured, simulating the joint of a surgical instrument. The removal of all test soil indicates effective cleaning.



THE ONLY RELIABLE METHOD TO ENSURE THE WASHER-DISINFECTOR IS WORKING ADEQUATELY.



Clear Instrument Tray System



Clear Instrument Trays		1 EACH	2 EACH	3+ EACH
CS-T1.6	Small 4" x 6½" x 7" (1.6 qt)	\$69.00/ea	\$64.00/ea	\$59.00/ea
CS-T6.3	Large 4" x 11" x 13" (6.3 qt)	\$109.00/ea	\$105.00/ea	\$99.00/ea
CS-T20.6	Extra Large 6" x 12¾" x 21" (20.6 qt)	\$179.00/ea	\$174.00/ea	\$169.00/ea

These instrument trays are clear, completely see-thru hard non-toxic polycarbonate material that include strainers and covers. Use for soaking and/or storage and transporting of instruments and devices. Special labels may be used to identify the contents and are easy peel-off for changing expiration dates and other information.

Self-Straining Instrument Soaking Trays



Self-Straining Instrument Soaking Trays		1 EACH	2 EACH	3+ EACH
CS-TH10	Large 10" x 5" x 2½" (LxWxD)	\$129.00/ea	\$124.00/ea	\$119.00/ea

The strainer lifts automatically when lid is opened. Each tray comes with it's own counter, which can be reset for each change of HLD solution! Recommended for use with glutaraldehyde solutions (not compatible with phenolic solutions). Available in 2 sizes.

Sterilizers & Accessories To Meet Your Needs.



Sterilizer Door
Gaskets
page 24

Sterilizer
Cleaners
page 25

Water
Distillers
page 26

Refurbished
Sterilizers
page 27

STERILIZER DOOR GASKETS



Contact a representative today to determine which gasket is compatible with your sterilizer. Verify with your sterilizer manufacturer how often they recommend replacing the door gasket on your sterilizer model. It is one essential and critical part of ensuring sterility of your instruments.

Don't see your sterilizer model listed here? Call 800.701.0186 to see if your door gasket is in stock!

DID YOU KNOW

Leaky door gaskets are a common cause of failed spore tests?

STERILIZER TIP: We recommend changing your sterilizer door gasket minimally of once a year to ensure there is no pressure loss in the chamber, which causes sterilizer failure. Check your sterilizer manual for manufacturer's recommendations.

Door Gaskets

	CROSS-SECTION	STOCK NO.	STERILIZER MAKE & MODEL	OEM#	SIZE	1 EACH
MIDMARK		DG-M7	SpeedClave® 7 & M7	H98134, 33166 & 002-0243-00	7.5" OD x .22" diameter C/S (19.1cm x .55cm OD)	\$54.00/ea
		DG-M9	M9, M9D	053-0366-00	~10.75" OD (~27.3cm OD)	\$89.00/ea
		DG-M11	M11 & M11D	053-0527-00	~12.86" OD (~32.7cm OD)	\$99.00/ea
SCICAN		DG-S2000	Statim 2000 - replacement cassette seal	01-100028S	~7.75" x 11.75" (~19.7cm x 29.9cm)	\$109.00/ea
		DG-S5000	Statim 5000 - replacement cassette seal	01-101649S	~7.75" x 15.75" (~19.7cm x 40.0cm)	\$119.00/ea
PELTON & CRANE		DG-OCM	OCM®	004014	8.25" OD (21.0cm OD)	\$79.00/ea
		DG-OCR	OCR® & OCR+®	004341	10.25" OD (26.1cm OD)	\$85.00/ea
		DG-PC8	Delta 8" (AE), Delta Q8, Delta XL8, Validator 8, Validator Plus 8	1539241 & 019613	9.88" OD (25.1cm OD)	\$69.00/ea
		DG-PC10	Delta AE8, XL8, Q8, AF10, XL10 & Q10, Validator Plus 10	1539407 & 019763	11.31" OD (28.7cm OD)	\$85.00/ea
PRESTIGE		DG-KV	KavoClave & Prestige 2100 Series Old Verison I, New Version II, & Prestige 2100 Classic Version III	219500	9.66" OD (24.5cm OD)	\$129.00/ea
TUTTNAUER		DG-T1730	Valueklave 1730, 1730E/EK/M/MK	02610020	~7.37" OD (~18.7cm OD)	\$69.00/ea
		DG-T2340	EZ9 & 2340 E/EA/EK/EKA/M* (*Serial No. 8806 and higher)	02610118	~9.50" OD (~24.1cm OD)	\$99.00/ea
		DG-T2340M	2340M (for Serial No. 8805 and below)	02610005	~9.422" OD (~23.9cm OD)	\$89.00/ea
		DG-T2540	EZ10 & 2540 E/EA/EK/EKA/M/MK	02610023	~10.59" OD (~26.9cm OD)	\$99.00/ea
		DG-T3870	3870E/EA/EHS/M	02610019	~15.95" OD (~40.5cm OD)	\$169.00/ea

Duraclean Sterilizer Cleaner

Duraclean Autoclave Cleaner		1 BOTTLE	2 BOTTLES	3+ BOTTLES
AC-DC1	32 Ounce Bottle, 1 Bottle	\$35.00/btl	\$32.00/btl	\$29.00/btl

32 oz. bottle concentrate. Used to clean grease and grime from autoclaves and retards the accumulation of scale deposits. Mix 6 oz. with 64 oz. of distilled water. One 32 oz. bottle per 5 cleanings. Compares to Omni-Cleaner™.



Chamber Brite™ Sterilizer Cleaner

Chamber Brite™ Sterilizer Cleaner		1 BOX	2 BOXES	3+ BOXES
AC-CB	Chamber Brite™, 10 packs/box	\$95.00/bx	\$92.00/bx	\$89.00/bx

Manufactured by Tuttnauer™, these powdered autoclave cleaner come in pre-measured packs. They are designed specifically for Gravity Steam sterilizers. Effectively cleans and de-scales the reservoir chamber. Easy-to-use! Chamber Brite™ is not intended for use with any pre/post vacuum sterilizers.



Tutt-Clean™ Sterilizer Cleaner

TUTT-CLEAN™ Sterilizer Cleaner		1 BOX	2 BOXES	3+ BOXES
AC-TC	TUTT-CLEAN™ Cleaner, 10 packs/bx	\$55.00/bx	\$52.00/bx	\$49.00/bx

Designed for use in gravity-type autoclaves. TUTT-CLEAN™ is a mild acidic salt cleansing compound formulated for descaling and removal of water deposits, oxides and other sediments from sterilizers. Always follow equipment manufacturer's instructions and recommended maintenance schedules. TUTT-CLEAN™ is not intended for use with any pre/post vacuum sterilizers.

Note: Comparable to Chamber Brite™, by Tuttnauer™.



Mid-Clean™ Sterilizer Cleaner

Mid-Clean™ Sterilizer Cleaner		1 BOTTLE	2 BOTTLES	3+ BOTTLES
AC-MC	Mid-Clean™ Cleaner, 16 oz. bottle	\$39.00/btl	\$36.00/btl	\$33.00/btl

MID-CLEAN™ is a multi-purpose, slightly alkaline, mild solution developed to remove grease and grime from sterilizers. Used as directed, MID-CLEAN™ will also retard the accumulation of scale deposits. Always follow equipment manufacturer's instructions and recommended maintenance schedules. MID-CLEAN™ is intended to be used in all Midmark® sterilizers.



Sci-Dry™ Sterilizer Cleaner

SCI-DRY™ Statim Cleaner		1 BOTTLE	2 BOTTLES	3+ BOTTLES
AC-SD2	Sci-Dry™ 2 oz. spray bottle	\$40.00/btl	\$37.00/btl	\$34.00/btl
AC-SD16	Sci-Dry™ 16 oz. refill bottle	\$71.00/btl	\$68.00/btl	\$65.00/btl

SCI-DRY™ is designed to protect stainless steel instruments and will promote even coating of water on the internal surfaces of the cassette without beading. This will reduce spotting on the cassette as well as on the instruments. SCI-DRY™ is intended to be used with all SciCan® cassette sterilizers including models 2000, 5000 and 7000.



Premium Water Distiller



Premium Water Distiller		1 EACH
WD-80	Premium Water Distiller, 3 litres, 15 yr warranty	\$1,095.00/ea
WD-80-C	Additional water storage container, stores 3 litres	\$39.00/ea


The solution for all small to mid-sized offices or facilities requiring less than 3 gallons of distilled water per day. It is easy-to-use (operates much like a coffee maker) and there is no installation required. Simply fill up the boiling tank with tap water, turn the appliance on, and it will produce approximately 3 litres of distilled water per cycle! No filters to replace, and the only maintenance is a vinegar cycle once per month. The lifespan of this distiller is significantly longer than other manual distillers making it an excellent investment that will pay you back many times over.

KEY FEATURES:

- Eliminates water contaminants
- Protects your equipment
- 800 watts of power, plugs into any standard 120v outlet
- Produces 3L of distilled water per cycle
- Stainless steel
- Timer shut-off
- Built-to-Last: 15 year limited warranty
- Made in the U.S.A.

SPECS:

- Width 13.75" (35 cm)
- Depth 9.875" (25 cm)
- Height 16.5" (42 cm)
- Machine Weight 26 lbs. (12 kg)
- Current Draw 3.5 amps
- Total Wattage 800 watts



*On all stainless steel components



Using bottled water?

Or a low-quality distiller?

If your office uses either of these methods, you are most likely spending too much money. Our Premium distiller offers the most cost-effective solution for distilled water. See comparison below...

VALUE COMPARISON*

WD-80 DISTILLER

15 year cost: \$895
15+ year lifespan
\$895 x 1 (15 yr lifespan)
= \$895



OTHER UNITS

15 year cost: \$3,600
Approx. 2 year lifespan
\$480 x 7.5 (2 yr lifespan)
= \$3,600



BOTTLED WATER

15 year cost: \$6,300
Ongoing
\$35/month x 15 yrs
= \$6,300



* Enter in your own costs to calculate actual savings!

**6 MONTH
WARRANTY**

on all refurbished
sterilizers and
autoclaves

*Need To Replace Your Current Sterilizer(s) But Don't
Want To Pay The Price for A Brand New Unit?*

**Refurbished Sterilizers offer the same quality and performance
from the brands you trust, at a better savings for you.**

REFURBISHED STERILIZERS

Refurbished Sterilizers

Contact us for refurbished Steam and Dry Heat Sterilizers!
All used sterilizers are refurbished by a qualified technician
- for up to half the price of a brand new unit.

CALL US FOR SPECIAL PRICING ON

- Pelton & Crane
- Tuttnauer
- Cox
- Statim
- Midmark
- And other leading brands!

What does a Refurbished Sterilizer include?

- Complete Disassembling of the Autoclave
- Defective and/or Questionable Parts disposed of
- Frame/Chassis sanded and painted if needed
- Chamber cleaned free of dirt & mineral buildup
- Reservoir cleaned to perform like new
- All O-Rings, Gaskets, Seals, Filters replaced
- Valves & Switches replaced as needed
- Complete Re-assembling
- Hinges, Latches, moving parts are lubricated
- Wire checked for shorts or ground leakage
- Plumbing inspected for water or pressure leaks
- Put thru QC (Quality Control)
 - Tested by Factory Trained Technician
 - Calibrated to Mfg. Specifications
 - Closely inspected for any abnormalities
 - Spore Tested
 - Labeled and Documented
- Fitted with Rack, Trays, & Operators Manual
- Bubble wrapped & double boxed for shipment



Call 800.701.0186 for more information! Have a question or need more options?
Email info@advantechsterilizers.ca and our sales team will be happy to help!

- www.advantechsterilizers.ca -



Toll Free 800.701.0186
info@advantechsterilizers.ca
www.AdvanTechSterilizers.ca