



CLAVE45+
AUTOCLAVE
STEAM STERILIZER 45L - CLASS S



LARGE CAPACITY & EFFICIENT FOR ALL APPLICATIONS


- Easy to Repair – Quality Parts Made from Europe
- Built in Printer Included
- Double Door Protection System – Prevents door from opening during cycle
- Closed Door Drying
- Automatic Door Closing Mechanism
- Dry-Assist Technology for Fast Drying with post-vacuum
- Perfect for Labs, Dental, Medical, Veterinary & Food Processing Facilities



Trolley Included with Each Unit

| STERILIZATION PROGRAMS | | | |
|------------------------|---------|---------------|--|
| Programs | Temp °C | Pressure k.pa | Type |
| UNWRAPPED | 134 | 210 | Unwrapped solid material |
| | 121 | 110 | |
| WRAPPED | 134 | 210 | Wrapped solid or hollow material |
| | 121 | 110 | |
| LIQUID | 121 | 110 | Liquid |
| PRION | 134 | 210 | Unwrapped porous material |
| | | | Single-wrapped porous material |
| | | | Dual-wrapped porous material |
| | | | Single-wrapped hollow material |
| | | | Dual-wrapped solid and hollow material |

| FLIGHT CLAVE45+ | |
|-----------------------|-----------------------------------|
| Capacity | 45L |
| Load Size | 22lbs-Wrapped |
| Chamber Size | Ø12.5"x24.3"(45L) |
| Cycle Time | 70min Wrapped |
| Voltage | 220V - 15A Breaker 2600W Power |
| Internal Memory | 10,000 Cycles |
| Printer | Built in (Included) |
| Datalogger | USB |
| Warranty | 2 Year Parts Warranty |
| Trays | 2 or 5 |
| Clean Water Reservoir | Internal Tank |
| Waste Water | Internal Tank and Auto Drain |
| Bowie Dick Test | Not Required |
| Class of Sterilizer | Class S |



| | |
|---------------------|-------|
| W | 25.2" |
| H | 22" |
| D | 25.2" |
| Door Open Width(WO) | 29.9" |
| Door Open Depth(HO) | 49.8" |