




12
MINUTE
CYCLE

CLAVE8+
AUTOCLAVE
ULTRA FAST CLASS S






SPECIFICATIONS

-  12 Min Cycle
-  Small Footprint
-  Big Power





**Standard Configuration
2 Trays (3 trays Optional)

USB Datalogger
& Internal Printer Included



Model	Clave8+
Voltage	110V
Rated Power	1750W
Chamber Size	Ø6.7"x12.6"(8L)
Weight	82 lbs

DIMENSIONS	
W	16.5"
H	14.6"
D	20.7"
L1	14"
L2	13"
Td	10.4"
Tw	5.1"



Waste Water Drain Distilled Water Drain