

FEATURES

- 2 High speed Connections
- 1 Low speed Connection
- Easy Fill Oil Reservoir
- Clean and mess-free oil sponge
- Compressed Air Input with Easy connection
- Individual programmed Cycles that provides thorough cleaning and lubrication
- Easily detachable door for easy cleaning

PROGRAMMED CYCLES

S: Short Cycle for High Speed Handpieces
(Lubricating and Cleaning Time: 35 seconds)

L: Long Cycle for Low Speed Handpieces
(Lubricating and Cleaning Time: 45 seconds)

E: Extra Long Cycle for Speed Reduction
Handpieces (Lubricating and
Cleaning Time: 50 seconds)



EASY TO CLEAN

The front door can be completely removed to allow fast, easy cleaning of the system



OIL RESERVOIR

The oil reservoir has a 500ml capacity, reducing the need for frequent refilling. The clean easy to use fill tank also has an oil level indicator to allow the user to see when to add more oil.

ACCOMMODATES ALL BRANDS OF HANDPIECES

All brands of E-Type low speed handpieces connect directly to the coupling supplied. The highspeed handpiece connects directly to a 4-hole connection or a quick disconnect coupler can be installed.



Technical Specifications

Rated Voltage	AC120V +/-10% 60Hz
Input Voltage	120V: 19VA
Air Pressure	0.35-0.6MPa - (50-85psi)
Dimensions	W 280mm x D275mm x H 360mm W 11" x D10.8" x H 14.2"
Weight	6.8kg or 15lbs